

ECONOMICAL, ONLY
DISPOSE OF THE LINER!

Disposable
Polypropylene
Liners




Chill Trough™

Reagent Reservoir

With disposable polypropylene liners for multi-channel pipettors. ChillTrough™ is designed to keep reagents cold at 32 - 40°F during sample preparation. Supports a disposable polypropylene trough liner. Coolant is sealed within the trough cavity and is environmentally safe.

- Easy to use:
 - Store the ChillTrough™ in a freezer.
 - Remove from the freezer before use.
 - Drop-in disposable polypropylene liner into the ChillTrough™.
- Save money:
 - Liners are made of thin, disposable lightweight polypropylene.
 - Troughs can be reused.

COOLING REAGENT RESERVOIRS



**Internal Coolant
Keeps Reagents
Cold During
Preparation!**

Part #	Description	Items/pk
CTKIT-100	Reagent Reservoir Kit for multi-channel pipettors: Contains 2 regular troughs, 2 ChillTroughs™, 400 disposable polypropylene trough liners, 100 mL capacity	1 Kit
CT-L-100	Disposable polypropylene ChillTrough™ liners , 100 mL capacity	400
CT-RG-100	Regular Trough , polypropylene, 100 mL capacity	6
CT-CT-100	ChillTrough™ polypropylene reagent reservoirs , 100 mL capacity	6