



Polypropylene

**Anti-Static
Polypropylene**

**Borosilicate
Glass**

Disposable Weighing Funnels Patent No. 6179022 Patent No. 415662

- Sit flat on balance - No Rolling, No Spilling!
- Pouring necks designed to fit volumetric flasks, micro-tubes or any narrow neck apparatus for spill free transfer
- Residual sample can be rinsed to guarantee 100% transfer
- Suitable for safe handling of hazardous materials
- Available in 5 different sizes

Part #	Weighing Funnel Capacity (mL)	Fits Volumetric Flask Volume (mL)	Items/pk
--------	----------------------------------	--------------------------------------	----------

DISPOSABLE ANTI-STATIC POLYPROPYLENE

ASWFXS	2 mL solids/1 mL liquid	96 Well Microtubes	96
ASWF1S	2 mL solids/2 mL liquid	5, 10, 25, 50	50
ASWF2M	10 mL solids/7 mL liquid	100, 200, 250	50
ASWF3L	20 mL solids/10 mL liquid	500, 1000	20
ASWF4XL	60 mL solids/25 mL liquid	1000, 2000	10

DISPOSABLE POLYPROPYLENE

DPWF-PP-XS	2 mL solids/1 mL liquid	96 Well Microtubes	96
DPWF-PP1-S	2 mL solids/2 mL liquid	5, 10, 25, 50	50
DPWF-PP2-M	10 mL solids/7 mL liquid	100, 200, 250	50
DPWF-PP3-L	20 mL solids/10 mL liquid	500, 1000	20
DPWF-PP4-XL	60 mL solids/25 mL liquid	1000, 2000	10

BOROSILICATE GLASS

GWF-BS-XS	2 mL solids/1 mL liquid	96 Well Microtubes	6
GWF-BS1-S	2 mL solids/2 mL liquid	5, 10, 25, 50	6
GWF-BS2-M	10 mL solids/7 mL liquid	100, 200, 250	6
GWF-BS3-L	20 mL solids/10 mL liquid	500, 1000	6
GWF-BS4-XL	60 mL solids/25 mL liquid	1000, 2000	6