



TWD Scientific, LLC  
 10555 86<sup>th</sup> Avenue  
 Pleasant Prairie, WI 53158

888.323.3585 TEL | 262.605.3262 FAX

www.twdtradewinds.com

<b>PART #:</b>	<b>LOT #:</b>
<b>DF8988P-S</b>	<b>201906-01</b>
<b>DESCRIPTION:</b>	
<b>DISPOSABLE REINFORCED NYLON STERILE WHITE FORCEPS WITH BLUNT TIP.          FORCEPS ARE INDIVIDUALLY PEEL-PACKED.          *INTENDED FOR SINGLE USE ONLY*</b>	
<b>SPECIFICATIONS:</b>	<b>MATERIAL:</b>
<b>LENGTH: 122.5mm          WIDTH: 22.5mm          TIP OPENING: 15mm</b>	<b>REINFORCED NYLON</b>
<b>EXPIRATION DATE (MONTH/DAY/YEAR): 05/24/2024</b>	

## CERTIFICATE OF STERILIZATION

The validation and routine control of the sterilization process are carried out in accordance with the following standards:

ISO 11135: Sterilization of health care products – Ethylene oxide – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices

Provided appropriately stored, products remain sterile for at least 5 years from the sterilization date.

## CERTIFICATE OF CONFORMITY

We hereby certify that product mention above meets the specifications stated.

Product: Sterile Forceps

Lot #: 201906-01

Expiration (Month/Day/Year): 05/24/2024