



TWD Scientific, LLC  
10555 86<sup>th</sup> Avenue  
Pleasant Prairie, WI 53158

888.323.3585 TEL | 262.605.3262 FAX

[www.twdtradewinds.com](http://www.twdtradewinds.com)

PART #:	LOT #:
DF8988P-SPT	201906-02
DESCRIPTION:	
DISPOSABLE REINFORCED NYLON STERILE WHITE FORCEPS WITH POINTED TIP. FORCEPS ARE INDIVIDUALLY PEEL-PACKED. <small>*INTENDED FOR SINGLE USE ONLY*</small>	
SPECIFICATIONS:	MATERIAL:
LENGTH: 110mm WIDTH: 32.5mm TIP OPENING: 25mm	REINFORCED NYLON
EXPIRATION DATE (MONTH/DAY/YEAR): 05/24/2024	

## CERTIFICATE OF STERILIZATION

The validation and routine control of the sterilization process are carried out in accordance with the following standards:

ISO 11135: Sterilization of health care products – Ethylene oxide – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices

Provided appropriately stored, products remain sterile for at least 5 years from the sterilization date.

## CERTIFICATE OF CONFORMITY

We hereby certify that product mention above meets the specifications stated.

Product: Sterile Forceps

Lot #: 201906-02

Expiration (Month/Day/Year): 05/24/2024