

CERTIFICATE OF CONFORMITY

PART #:	LOT #:
DPWF-PP2-M	20231005
DESCRIPTION:	
MEDIUM DISPOSABLE POLYPROPYLENE WEIGHING FUNNEL. FITS INTO 100-250ML VOLUMETRIC FLASKS. <small>*INTENDED FOR SINGLE USE ONLY*</small>	
SPECIFICATIONS:	MATERIAL:
LENGTH: 62.5mm FUNNEL O.D.: 35mm POURING NECK O.D.: 10mm CAPACITY: 10mL SOLID/7mL LIQUID	VIRGIN POLYPROPYLENE

TWD Scientific, LLC certifies that:

1. The product will not be altered or modified in any manner from the original condition after the product(s) are manufactured and packaged.
2. The product(s) meet the description specifications as detailed in the product description above.

Prepared by: Christina Weber, Customer Service

