

CERTIFICATE OF CONFORMITY

PART #:	LOT #:
DPWF-PP1-S	20240404
DESCRIPTION:	
SMALL DISPOSABLE POLYPROPYLENE WEIGHING FUNNEL. FITS INTO 5-50ML VOLUMETRIC FLASKS. <small>*INTENDED FOR SINGLE USE ONLY*</small>	
SPECIFICATIONS:	MATERIAL:
LENGTH: 43.2mm FUNNEL O.D.: 22.9mm POURING NECK O.D.: 7.6mm CAPACITY: 2mL SOLID/2mL LIQUID	VIRGIN POLYPROPYLENE

TWD Scientific, LLC certifies that:

1. The product will not be altered or modified in any manner from the original condition after the product(s) are manufactured and packaged.
2. The product(s) meet the description specifications as detailed in the product description above.

Prepared by: Christina Weber / Customer Service

