

CERTIFICATE OF CONFORMITY

PART #:	LOT #:
ASWFXS	20220503A-12A
DESCRIPTION:	
EXTRA SMALL DISPOSABLE ANTI-STATIC POLYPROPYLENE WEIGING FUNNEL. FITS INTO MICRO-TUBES AND 96 DEEP WELL PLATES. <small>*INTENDED FOR SINGLE USE ONLY*</small>	
SPECIFICATIONS:	MATERIAL:
LENGTH: 37.5mm FUNNEL O.D.: 15mm POURING NECK O.D.: 5mm CAPACITY: 2mL SOLID/1mL LIQUID	VIRGIN POLYPROPYLENE

TWD Scientific, LLC certifies that:

1. The product will not be altered or modified in any manner from the original condition after the product(s) are manufactured and packaged.
2. The product(s) meet the description specifications as detailed in the product description above.

Prepared by: Christina Weber / Customer Service

