

# CERTIFICATE OF CONFORMITY

<b>PART #:</b>	<b>LOT #:</b>
DPWF-PP3-L	20220506
<b>DESCRIPTION:</b>	
LARGE DISPOSABLE POLYPROPYLENE WEIGHING FUNNEL. FITS INTO 500-1000ML VOLUMETRIC FLASKS <small>*INTENDED FOR SINGLE USE ONLY*</small>	
<b>SPECIFICATIONS:</b>	<b>MATERIAL:</b>
LENGTH: 88.90mm FUNNEL I.D.: 48.26mm POURING NECK O.D.: 17.78mm CAPACITY: 20mL SOLID/10mL LIQUID	VIRGIN POLYPROPYLENE

TWD Scientific, LLC certifies that:

1. The product will not be altered or modified in any manner from the original condition after the product(s) are manufactured and packaged.
2. The product(s) meet the description specifications as detailed in the product description above.

Prepared by: Christina Weber / Customer Service

 10/27/23